

technoclon

TECHNOSCREEN® ADAMTS-13 Activity

Semiquantitative flow through assay for the determination of ADAMTS-13 activity levels






Your partner for ADAMTS-13

PRODUCT INFORMATION

The TECHNOSCREEN® ADAMTS-13 is a semiquantitative flow through assay for the determination of ADAMTS-13 activity levels in human citrated plasma. This assay is intended as a screening tool for estimating ADAMTS-13 activity.

Patient plasma is incubated with a vWF fragment. ADAMTS-13 specific cleavage of the vWF is then detected with a cleavage site specific antibody by the amount of color development. The amount of cleavage directly correlates to the ADAMTS-13 activity level in the plasma sample.



CORRELATION

200 patient samples with different ADAMTS-13 activity levels were tested both using the TECHNOSCREEN® ADAMTS-13 Activity Test and the TECHNOZYM® ADAMTS-13 Activity ELISA .


Visual interpretation was performed by multiple independent investigators.

Based on a cut off of 0.1 IU/ml:


positive predictive value	99%
negative predictive value	80%
sensitivity	94%
specificity	97%


TECHNOSCREEN® ADAMTS-13

EASY & FAST




Short instruction with a pre-incubation tube holder included in every package





Large and small pipettes for the pre-incubation mixture included

Easy test procedure for rapid screening with ready to use reagents



For the interpretation of the results a color chart is included in the kit

No additional instrument required

Reduce time to diagnose

	Product	Description	REF	Package
K	TECHNOSCREEN® ADAMTS-13 Activity	10 x test devices 10 x pre-incubation tubes 1 x 4 ml R1 1 x 1 ml R2 1 x 1 ml R3 11 x large pipettes (300 µl) 22 x small pipettes (60 µl) 1 x Short instruction pad 1 x Reference color chart	5700100	10 tests
	Demonstration video for TECHN	NOSCREEN® ADAMTS-13 Activity		



available on the TECHNOCLONE Youtube Channel

QUANTITATIVE KITS AND KITS FOR FURTHER ADAMTS-13 DIAGNOSIS

Product	Description	REF	Package
TECHNOZYM® ADAMTS-13 Activity ELISA Chromogenic ELISA	12 x 8 anti GST strips 1 x 12 ml HRP Conjugate 1 x 30 ml Reaction Buffer 1 x 53 ml Wash Buffer Concentrate 1 x 12 ml TMB Substrate 2 x 6 ml GST-vWF73 Substrate 1 x 12 ml Stop Solution 8 x 0.5 ml Calibrators and Controls	5450701	96 tests
TECHNOZYM® ADAMTS-13 Antigen ELISA Chromogenic ELISA	12 x 8 anti ADAMTS-13 strips 1 x 0.3 ml anti-ADAMTS-13 Conjugate 1 x 80 ml Wash Buffer Concentrate 1 x 90 ml Incubation Buffer 1 x 12 ml TMB Substrate 1 x 12 ml Stop Solution 7 x 0.5 ml Calibrators and Controls	5450601	96 tests
TECHNOZYM® ADAMTS-13 INH ELISA Chromogenic ELISA	12 x 8 rec. ADAMTS-13 strips 1 x 0.3 ml anti-IgG Conjugate 1 x 80 ml Wash Buffer Concentrate 1 x 90 ml Incubation Buffer 1 x 12 ml TMB Substrate 1 x 12 ml Stop Solution 7 x 0.5 ml Calibrators and Controls	5450401 5450451	96 tests 48 tests
TECHNOZYM® ADAMTS-13 Activity/Antigen ELISA Fluorogenic ELISA	12/6 x 8 anti-ADAMTS-13 strips 1 x 0.3 ml anti-ADAMTS-13 Conjugate 1 x 80 ml Wash Buffer Concentrate 1 x 90 ml Incubation Buffer 2/1 x 3 ml Activity Substrate 1 x 6 ml Antigen Substrate 1 x 0.7 ml Stable Peroxidase Solution 1 x 6 ml Stop Solution 7 x 0.5 ml Calibrators and Controls	5450501 5450551	2 x 48 tests 48 tests

Technoclone Herstellung von Diagnostika und Arzneimitteln GmbH

Brunner Str. 67 | 1230 Vienna | Austria

Phone: +43 1 86373-0 Fax: +43 1 86373-44

sales@technoclone.com www.technoclone.com















