

2016/2017

Hospital Trolleys



SST Trolleys

page 1-8

SST Accessories

pages 9-10

GFS Trolleys

pages 11-18

GFS Accessories

page 19

ALS Trolleys

pages 20-26

ALS Accessories

pages 27-30



BARWORTH
MEDICAL 

Barworth Medical (UK) Ltd, Trackat House, 24 Turnpike Road
Blunsdon, Swindon, Wilts SN26 7EA

Tel: 01793 729229 Fax: 01793 728820

email: info@barworthmedical.co.uk

www.barworthmedical.co.uk

SST Trolleys

10525

Hospital Trolley

- Stainless steel structure
- Lower stainless steel shelf, with rails on three sides
- Upper ABS platform on top with extended up rims on three sides
- Four 125mm rolling castors, two lockable
- Bumper protection



L62 W48 H96 12

£430

exc VAT & delivery



10530

Hospital Trolley

- Stainless steel structure
- Middle and lower stainless steel shelves, with rails on three sides
- Upper ABS platform on top with extended up rims on three sides
- Four 125mm rolling castors, two lockable
- Bumper protection



L62 W48 H96 16

£495

exc VAT & delivery



10810

Hospital Trolley

- Stainless steel structure
- Middle and lower stainless steel shelves, with rails on three sides
- Upper ABS platform on top with extended up rims on three sides
- Four 125mm rolling castors, two lockable
- Bumper protection



L100 W48 H96 22

£730

exc VAT & delivery



10800

Hospital Trolley

- Stainless steel structure
- Lower stainless steel shelf, with rails on three sides
- Upper ABS platform on top with extended up rims on three sides
- Four 125mm rolling castors, two lockable
- Bumper protection



L100 W48 H96 18

£595

exc VAT & delivery







SST Trolleys

14170

Single width Trolley



- Steel chassis, durable white finish
- PVC coated upper shelf with up-rim on three sides
- PVC coated lower shelf
- Single drawer
- 80mm wheels

			
L48	W47	H85	14

£330
exc VAT & delivery

14150

Double width Trolley



- Steel chassis, durable white finish
- PVC coated upper shelf with handrail on three sides
- PVC coated lower shelf
- Single drawer
- 80mm wheels





			
L62	W48	H85	16

£330
exc VAT & delivery

10510

Hospital Trolley

- Stainless steel structure
- Upper ABS platform on top with extended up rims on three sides
- PVC drawers with automatic shutting system
- Lower stainless steel shelf with rails on three sides
- Four 125mm rolling castors, two lockable
- Bumper protection

			
L62	W48	H96	18





£746
exc VAT & delivery



10555

Hospital Trolley






- Stainless steel structure
- Middle and lower stainless steel shelves, both with rails on three sides
- Upper ABS platform on top with extended up rims on three sides
- PVC drawer with automatic shutting system
- Four 125mm rolling castors, two lockable
- Bumper protection

			
L62	W48	H96	23

£665
exc VAT & delivery



Colour Options

	
Blue	Lt Green
	
Purple	Red
	
White	Orange

SST Trolleys

10610

Multifunctional Hospital Trolleys

- Stainless steel structure
- Upper ABS platform on top with extended up rims on three sides
- Drawers with automatic shutting system, 110% total opening
- Lower open compartment
- 125mm rolling wheels, two with brakes
- Identification card holder
- Bumper protection



10535	62	53	100	33	0,33
10540	62	53	100	30	0,33
10605	62	53	100	36	0,33
10610	62	53	100	40	0,33

£1175

exc VAT & delivery



10610

10535
£980

exc VAT & delivery

10540

£1055

exc VAT & delivery



10540

10605

£1125

exc VAT & delivery



10605



10535

Colour Options



SST Trolleys



10550

10550/10565

Multifunctional Hospital Trolleys

- Stainless steel structure
- Upper ABS platform on top with extended up rims on three sides
- Drawers with automatic shutting system, 110% total opening
- 125mm rolling wheels, two with brakes
- Identification card holder
- Bumper protection



10550	62	53	100	47	0,33
10565	62	53	100	32	0,33

10550

£1549

exc VAT & delivery

10565


£1035

exc VAT & delivery



10565

Colour Options








SST Trolleys

10615/10500

Multifunctional Hospital Trolleys

- Stainless steel structure
- Upper ABS platform on top with extended up rims on three sides
- Drawers with automatic shutting system, 110% total opening
- 125mm rolling wheels, two with brakes
- Identification card holder
- Bumper protection

					
10500	62	53	100	40	0,33
10615	62	53	100	44	0,33

10615

£1329

exc VAT & delivery


10500

£1240

exc VAT & delivery



Colour Options



SST Trolleys



10760

Multifunctional Hospital Trolleys

- Stainless steel structure
- Upper ABS platform on top with extended up rims on three sides
- Middle and lower stainless steel shelves both with guard rails on three sides
- Drawers with automatic shutting system, 110% total opening
- 125mm rolling wheels, two with brakes
- Identification card holder
- Bumper protection



	100	48	100	44	0,48
10750	100	48	100	44	0,48
10760	100	48	100	48	0,48
10780	100	48	100	42	0,48

£1345

exc VAT & delivery

10750

£1210

exc VAT & delivery




10780

£1045

exc VAT & delivery

Colour Options



SST Trolleys

10575 + 10182/2 + 10182/13

Dressing and Treatment Trolleys

- Stainless steel structure
- Upper ABS platform on top with extended up rims on three sides
- Drawers with automatic shutting system, 110% total opening
- Identification card holder
- 125mm rolling wheels, two with brakes
- Bumper protection
- Equipped with accessories:
Ret.10182/2 Lateral tray retractable
Ret.10182/13 PVC waste bin



10575	62	53	100	38	0,33
10590	62	53	100	18	0,33

£1276

exc VAT & delivery



10575



10590

10590

£895

exc VAT & delivery

Colour Options



SST Trolleys

10850

Multifunctional Hospital Trolleys

- Epoxy steel coating structure
- Upper ABS platform on top with extended up rims on three sides
- Drawers with automatic shutting system, 110% total opening
- 125mm rolling wheels, two with brakes
- Identification card holder
- Bumper protection



10830	62	53	100	34	0,33
10840	62	53	100	36	0,33
10850	62	53	100	40	0,33

£1055

exc VAT & delivery



10850

10830

£995

exc VAT & delivery



10830

10840

£1115

exc VAT & delivery



10840

Colour Options



SST Trolley Accessories

10182-10

Oxygen Tank Holder

£68

exc VAT & delivery




10182-9

Defibrillator Monitor Shelf

£164

exc VAT & delivery



10182-6

IV Drip Rod Support, adjustable in height

£63

exc VAT & delivery



10182-25

Lock with Security Seal, including 100 security seals

£73

exc VAT & delivery



10182-11

Cardiac Back Board

£88

exc VAT & delivery




10182-70

Tilt out Bins Support

£172

exc VAT & delivery



10182-50

Bottles Support

£98

exc VAT & delivery



10182-3

Side Support

£50

exc VAT & delivery



10182-55

Gloves Box Support

£89

exc VAT & delivery



10182-28

Cable Guide Support

£165

exc VAT & delivery



10182-60

Gloves & Bottles Support

£128

exc VAT & delivery



SST Trolley Accessories


10182-4

Adjustable Divider for Drawer

- Flexible sliding divider with breaking points every 11mm

£55

exc VAT & delivery



10182-8

Base with multiple Sockets

£49

exc VAT & delivery



10182-30

Hand Sanitiser Dispenser Support

£177

exc VAT & delivery



10182-65

Laptop Support

£230

exc VAT & delivery



10182-26

Lateral Fixing Rail

£29

exc VAT & delivery




10182-1

Centralised Lock

£85

exc VAT & delivery



10182-13

PVC Waste Bin

£45

exc VAT & delivery



10182-2

Lateral Tray retractable

£76

exc VAT & delivery



10182-29

Sharp Disposal

£60

exc VAT & delivery



GFS Trolleys

70040/70045

Multifunctional Hospital Trolleys

- Aluminium/Steel epoxy coated frame
- ABS work surfaces with moulded up-rims on three sides
- Lateral push handle
- Both models have three equally spaced shelves
- 70040 has lower drawer with auto-close system
- Choice of drawer front/accent colour
- 100mm castor wheels, 2 braked



	65	56	94	26	0.34
70040	65	56	94	26	0.34
70045	65	56	94	18	0.34

70040

(with drawer)

£715

exc VAT & delivery

70045

(shelves only)

£570

exc VAT & delivery



70055/70060

Multifunctional Extra-Height Trolleys

- Additional height over standard trolleys
- Upper drawer + 2 shelves on 70055
- 4 equally spaced shelves on 70060



	65	56	128	27	0.46
70055	65	56	128	27	0.46
70060	65	56	128	20	0.46

70055

(with drawer)

£765

exc VAT & delivery

70060

(shelves only)

£695

exc VAT & delivery

Colour Options

drawer front,
side strips, wheels



Blue



Lt Green



Purple








Red

70000/70005L

Multifunctional Hospital Trolleys

- Aluminium/Steel epoxy coated frame
- ABS top surface with moulded up-rims on three sides
- Lateral push handle
- Upper drawer with auto-close system
- Choice of drawer front/accent colour
- 100mm castor wheels, 2 braked
- Low-level ABS shelf fitted to 70005L model

					
70000	65	56	94	24	0.34
70005L	65	56	94	24	0.34

70000

(no shelf)

£550

exc VAT & delivery

70005L

(with lower shelf)

£620

exc VAT & delivery



70000



70005L
Features a lower shelf to the same specifications as the upper surface



70035

70035

Multifunctional Hospital Trolley

- Aluminium/Steel epoxy coated frame
- ABS work surfaces with moulded up-rims on three sides
- Lateral push handle
- Three equally spaced ABS shelves
- Upper drawer with auto-close system
- Choice of drawer front/accent colour
- 100mm castor wheels, 2 braked

				
65	56	94	26	0.34

£695

exc VAT & delivery

Colour Options

drawer front,
side strips, wheels



Blue



Lt Green



Purple





Red

70005M

Multifunctional Hospital Trolley

- Aluminium/Steel epoxy coated frame
- ABS top surface, moulded up-rims on three sides
- Mid-position ABS shelf
- Lateral push handle
- Upper drawer with auto-close system
- Choice of drawer front/accent colour
- 100mm castor wheels, 2 braked

				
65	56	94	24	0.34

£620

exc VAT & delivery



70005M


Mid-Position
Shelf

GFS Trolleys



70050

Extra Height
Trolley



70050

Multifunctional Extra-Height Trolley

- Aluminium/Steel epoxy coated frame
- Additional height over standard trolleys
- ABS work surfaces with moulded up-rims on three sides
- Lateral push handle
- Slide-out writing shelf and upper drawer with auto-close system
- Choice of drawer front/accent colour
- 100mm castor wheels, 2 braked

65	56	128	28	0.46

£850

exc VAT & delivery

70020


Cure Trolley

- Aluminium/Steel epoxy coated frame
- ABS top surface with moulded up-rims on three sides
- Lateral push handle
- Three drawers with auto-close system
- Choice of drawer front/accent colour
- 100mm castor wheels, 2 braked

65	56	94	36	0.34

£725

exc VAT & delivery



70020
Shown with
Accessories

70240

Sanitizer
Dispenser Support

£177

exc VAT & delivery

70220

PVC
Waste Bin

£45

exc VAT & delivery

Colour Options

drawer front,
side strips, wheels



Blue



Lt Green



Purple








Red

70070/70075

Multifunctional Hospital Trolleys

- Aluminium and stainless steel structure
- Removable stainless steel upper and lower shelves
- High impact ABS base
- Lateral push handle
- Choice of single or twin drawer
- 80mm castor wheels, 2 braked
- Choice of drawer front/accent colour

					
70070	76	45	86	22	0.30
70075	76	45	86	22	0.30

70070

(single drawer)

£560

exc VAT & delivery

70075

(twin drawer)

£590

exc VAT & delivery



70250

Gauze Basins and Support
(Optional accessory)



Colour Options

drawer front,
side strips, wheels



Blue



Lt Green



Purple



Red








70065

Shown with
Accessories

70065

Multifunctional Hospital Trolley

- Aluminium and stainless steel structure
- Removable stainless steel upper and lower shelves
- High impact ABS base
- Full width drawer
- Fitted with PVC waste bin and basin holder
- Lateral push handle
- Choice of drawer front/accent colour
- 80mm castor wheels, 2 braked

				
76	45	86	22	0.30

70065

Trolley

£650

exc VAT & delivery

70250

Gauze basins + Support

£176

exc VAT & delivery






GFS Trolleys



70085/70090

Multifunctional Hospital Trolleys

- Aluminium and stainless steel structure
- High impact ABS base
- 70085 is wider, with two upper drawers, 3 removable stainless steel shelves, a PVC waste bin and lateral push handle
- 70090 has 2 removable stainless steel shelves
- 80mm castor wheels, 2 braked
- Choice of drawer front/accent colour

					
70085	90	45	94	43	0.38
70090	76	45	86	20	0.30

70085
£660
exc VAT & delivery






70090
£440
exc VAT & delivery



70115/70110

Multifunctional Hospital Trolleys

- Aluminium/Steel epoxy coated frame
- Wipe-clean laminate upper work surface
- Lower shelf is epoxy coated steel
- 70110 has full width upper drawer
- Choice of drawer front/accent colour
- 100mm castor wheels, 2 braked

					
70110	90	45	87	30	0.35
70115	90	45	87	21	0.35

70115
Trolley
£410
exc VAT & delivery

70110
Trolley with Drawer
£520
exc VAT & delivery

Colour Options

drawer front,
side strips, wheels



Blue



Lt Green



Purple



Red

70010/70015

Multifunctional Hospital Trolleys

- Aluminium and epoxy steel coating structure
- ABS upper top with extended up rims on three sides
- Lateral push handle
- High impact ABS base
- Drawers with automatic shutting system, 110% total opening
- 100mm rolling wheels, two with brakes
- Identification card holder



70010	65	56	94	30	0,34
70015	65	56	94	31	0,34

70010
£630

exc VAT & delivery

70015
£720

exc VAT & delivery



70010



70015
ABS Lower Shelf



70080

70080

Multifunctional Hospital Trolley with 3 Baskets

- Aluminium and epoxy steel coating structure
- ABS upper top with extended up rims on three sides
- Lateral push handle
- High impact ABS base
- Drawers with automatic shutting system, 110% total opening
- 100mm rolling wheels, two with brakes
- Identification card holder




70080	90	45	94	43	0.38
--------------	----	----	----	----	------

£920

exc VAT & delivery

Colour Options



GFS Trolleys

70030

Multifunctional Hospital Trolley

- Aluminium and epoxy steel coating structure
- ABS upper top with extended up rims on three sides
- Lateral push handle
- High impact ABS base
- Drawers with automatic shutting system, 110% total opening
- 100mm rolling wheels, two with brakes
- Identification card holder

					
70030	65	56	94	36	0,34

£895

exc VAT & delivery







70030

70025

Multifunctional Hospital Trolley

- Aluminium and epoxy steel coating structure
- ABS upper top with extended up rims on three sides
- Lateral push handle
- High impact ABS base
- Drawers with automatic shutting system, 110% total opening
- 100mm rolling wheels, two with brakes
- Identification card holder

					
70025	65	56	94	37	0,34

£820







exc VAT & delivery



70025

ABS Lower Shelf






Colour Options

	
Blue	Lt Green
	
Purple	Red
	
White	Orange

70100/70105

Multifunctional Hospital Trolleys

- Aluminium and stainless steel structure
- Stainless steel baskets
- High impact ABS base
- 100mm rolling wheels, two with brakes

					
70100	77	45	87	16	0.30
70105	77	45	87	16	0.30

70100

£660

exc VAT & delivery

70105


£610

exc VAT & delivery



70100

Three Grid Baskets 60 X 40 X 10cm



70105

Two Grid Baskets 60 X 40 X 20cm

Colour Options



Blue Lt Green



Purple Red



White Orange

GFS Trolley Accessories

70205

Cardiac Board

£80

exc VAT & delivery



70225

Sharp Disposal

£60

exc VAT & delivery



70220

PVC Waste Bin

£45

exc VAT & delivery



70200

Oxygen Tank Holder

110x780mm

£68

exc VAT & delivery



70200

Oxygen Tank Holder

140x570mm

£68

exc VAT & delivery




70245

Movable Separations

£45

exc VAT & delivery



70230

Lock with Security Seal,
including 100 security seals

£73

exc VAT & delivery



70215

Tilt out Bins
Support

£172

exc VAT & delivery




70210

IV Drip Rod Support,
adjustable in height

£63

exc VAT & delivery



70235

Laptop Support

£230

exc VAT & delivery



70240

Hand Sanitiser Dispenser Support

£177

exc VAT & delivery



ALS Trolleys



50005

Includes removable baskets
ISO modular 600 x 400mm



50005/50000

Multifunctional Hospital Trolleys

- Aluminium structure and stratified side panels with plasticised coating
- Upper ABS platform on top with extended up rims on three sides
- Base with bumper protection
- Drawers with automatic shutting system, 110% total opening
- Identification card holder
- Ergonomic pushing handle
- 125mm rolling wheels, two with brakes



50000	80	58	98	60	0.45
50005	80	58	98	65	0.45

50005

£1530

exc VAT & delivery


50000

£1373

exc VAT & delivery



Includes removable baskets
ISO modular 600 x 400mm



Colour Options



Blue Lt Green



Purple Red



White Orange



50000

ALS Trolleys

50030/50035

Multifunctional Hospital Trolleys

- Aluminium structure and stratified side panels with plasticised coating
- Upper ABS platform on top with extended up rims on three sides
- Base with bumper protection
- Drawers with automatic shutting system, 110% total opening
- Identification card holder
- Lower part with liner panels for baskets or shelves
- Ergonomic pushing handle
- 125mm rolling wheels, two with brakes



50030

Includes removable baskets
ISO modular 600 x 400mm



50030	80	58	98	60	0.45
50035	80	58	98	63	0.45

50030

£1125


exc VAT & delivery

50035

£1270


exc VAT & delivery

Includes removable baskets
ISO modular 600 x 400mm



50035

Colour Options



ALS Trolleys



50100

Includes removable baskets
ISO modular 600 x 400mm

50mm | 2

100mm | 6

50100/50110

Multifunctional Hospital Trolleys

- Aluminium structure and stratified side panels with plasticised coating
- Upper ABS platform on top with extended up rims on three sides
- Base with bumper protection
- Panel for baskets or shelves ISO 600 x 400mm
- Two doors
- Locking key
- Ergonomic pushing handle
- 125mm rolling wheels, two with brakes



50100	80	58	120	62	0.56
50110	80	58	156	75	0.72

50100

£1315

exc VAT & delivery

50110

£1555

exc VAT & delivery

Includes removable baskets
ISO modular 600 x 400mm

100mm | 8

200mm | 1

Colour Options



Blue Lt Green



Purple Red



White Orange



50110

50120/50115

Multifunctional Hospital Trolleys

- Steel structure with epoxy coating
- Stratified top
- Pushing handle
- Bumper protection
- 125mm rolling wheels, two with brakes



50115	74	42	110	25	0.34
50120	74	42	157	36	0.39

50120


£1075

exc VAT & delivery

50115


£730

exc VAT & delivery




50120

9 Guides for Baskets or Shelves
Includes removable baskets
ISO modular 600 x 400mm



50115

5 Guides for Baskets or Shelves
Includes removable baskets
ISO modular 600 x 400mm



ALS Trolleys



50230

50230

Multifunctional Hospital Trolley

- Aluminium structure and stratified side panels with plasticised coating
- Upper ABS platform on top with extended up rims on three sides
- Base with polyurethane integral bumper protection
- Two drawers
- Ergonomic pushing handle
- 125mm rolling wheels, two with brakes



50230	80	58	120	90	0.56
-------	----	----	-----	----	------

£1085

exc VAT & delivery

50220

Multifunctional Hospital Trolley

- One drawer for vertical clinical file folders
- 11 compartments for x-rays of 380 x 430mm
- Slide out writing desk



50220	80	58	120	77	0.56
-------	----	----	-----	----	------


£1165

exc VAT & delivery



50220

Colour Options



ALS Trolleys

50050/50055

Multifunctional Hospital Trolleys

- Aluminium structure and stratified side panels with plasticised coating
- Upper ABS platform on top with extended up rims on three sides
- Base with bumper protection
- Drawers with automatic shutting system, 110% total opening
- Identification card holder
- Lower part with liner panels for baskets or shelves
- Ergonomic pushing handle
- 125mm rolling wheels, two with brakes



50050



50050	80	58	98	70	0.45
50055	80	58	98	75	0.45

50050

£1240


exc VAT & delivery

50055


£1376

exc VAT & delivery

Includes removable baskets
ISO modular 600 x 400mm




Includes removable baskets
ISO modular 600 x 400mm



50055

Colour Options



Multi-Function Trolleys

H536

Multifunctional Hospital Trolley

- Upper and lower shelves with protective rim on three sides
- Two sliding drawers on metal glide runners with 110% opening
- Handles at either end
- Swivel castor wheels (100mm), 2 braked
- Main panels made of reinforced plastic polymer

85	43	84	19

*24kg with all accessories

£530

exc VAT & delivery



Optional Accessories

- Collection Bin:
Stainless steel 330 x 110 x 240mm
- End support with 3 containers
Stainless steel, container size
330 x 110 x 80mm
- Drawer Divider

Collection Bin

£76

exc VAT & delivery

End Support with 3 containers


£120

exc VAT & delivery

Drawer Divider

£60

exc VAT & delivery



▲ Optional drawer dividers

◀ Trolley fitted with all accessories

Drawer front colour options



Blue Green






Yellow Red

(Please state choice when ordering)

ALS Trolley Accessories

In addition to the trolleys, you will also find a set of trays and dividers available that are very easy to use, allowing you to personalise your compartments to your specific needs. You can compartmentalise the baskets both horizontally and transversely

- ISO Modular basket for drugs and expendable goods
- ISO Modular basket of ABS

			
50505	60	40	5
50510	60	40	10
50515	60	40	20
50520	40	30	10



50520
£28

exc VAT & delivery



50505
£30

exc VAT & delivery



50510
£35

exc VAT & delivery



50515
£40

exc VAT & delivery



ALS Trolley Accessories

Adjustable cross dividers and flexible sliding dividers with breaking points every 10mm




Adjustable cross dividers and flexible sliding dividers with breaking points every 10mm



Dividers for baskets of
600 x 400 x 50mm

50610




8 Compartment

£37

exc VAT & delivery

50605




15 Compartment

£56

exc VAT & delivery

50600



24 Compartment

£75

exc VAT & delivery

Dividers for baskets of
600 x 400 x 50mm

50615




12 Compartment

£33

exc VAT & delivery

50620




12 Compartment

£38

exc VAT & delivery

50625




12 Compartment

£44

exc VAT & delivery

Adjustable cross dividers and flexible sliding dividers with breaking points every 20mm



Adjustable cross dividers and flexible sliding dividers with breaking points every 10mm



Dividers for baskets of
600 x 400 x 100mm

50655




6 Compartment

£34

exc VAT & delivery

50650




8 Compartment

£43

exc VAT & delivery

50645




15 Compartment

£65

exc VAT & delivery

50660




9 Compartment

£44

exc VAT & delivery

50665




10 Compartment

£48

exc VAT & delivery

50670



12 Compartment

£58

exc VAT & delivery


ALS Trolley Accessories

50960

100 Security seals

£27

exc VAT & delivery



50900

Oxygen Tank Holder
110x780mm

£68

exc VAT & delivery



50970

Base with
multiple Socket

£49

exc VAT & delivery



50920

PVC Waste Bin

£45

exc VAT & delivery




50945

Side Basket of
400 x 300mm

£70

exc VAT & delivery




51010

Lateral Tray
retractable

£76

exc VAT & delivery




50940

Lateral pull-out Shelf

£115

exc VAT & delivery



50910

IV Drip Rod Support,
adjustable in height

£63

exc VAT & delivery



50955

Lock with Security Seal,
including 100 security seals

£73

exc VAT & delivery



ALS Trolley Accessories

50925

Sharp Disposal

£60

exc VAT & delivery



50905

Cardiac Back Board

£88

exc VAT & delivery



50975

Cable Guide Support

£165

exc VAT & delivery



50965

Lock with Key

£85

exc VAT & delivery




50935

Defibrillator
Monitor Shelf

£195

exc VAT & delivery



51005

Laptop Support

£230

exc VAT & delivery



50915

Tilt out Bins
Support

£172

exc VAT & delivery



50950

Side Support

£49

exc VAT & delivery





Barworth Medical (UK) Ltd, Trackat House, 24 Turnpike Road
Blunsdon, Swindon, Wilts SN26 7EA
Tel: 01793 729229 Fax: 01793 728820
email: info@barworthmedical.co.uk
www.barworthmedical.co.uk