

# CLINICON®

Laboratory Consumables & Equipment

## Introducing The New KapSafe™ Benchtop Recapper

Safely recap your sample tubes with a high-speed pneumatics free benchtop recapper. Recaps original vacuum tubes 13-16mm diameter, 75-115mm tall. KapSafe accommodates most analyser racks and provides multirack input for up to 20 racks.



### Technical data

Capacity:	Recaps up to 1,200 tubes per hour
Sample tubes:	Standard vacuum blood tubes
Caps:	Uses cap 2174-3056
Size:	60 x 60 x 83cms (Depth x Width x Height)
Voltage:	230VAC 50Hz / 110VAC 60Hz
Power:	60VA

✉ info@clinicon.co.uk

☎ +44 (0) 1273 463609

www.clinicon.co.uk