

manual pipette controllers

Profiller™ manual 435 / 436

Robust and easy to use, both models offer ergonomic shape and light weight. They provide smooth control in both aspirating and dispensing, using glass or plastic pipettes of volume up to 100 mL. One-year warranty.

Product advantages

- Large capacity vacuum bulb
- Soft action lever button for precise control
- Efficient blow-out system
- PTFE protection membrane filter
- Maintenance free
- Alternative colour coded nose cones



Squeezable large capacity bulbs



Precise control



Liquid blow-out



Protection filter



Colour cones

Ordering information - Instruments

| Description | Packaging | Cat. No. |
|--|-----------|----------|
| Profiller™ manual 435, egg style grey bulb | 1 / pk | 435.100 |
| Profiller™ manual 436, flat shaped blue bulb | 1 / pk | 436.100 |

Ordering information - Accessories

| Description | Packaging | Cat. No. |
|---|-----------|----------|
| PTFE Membrane filter, 0.45 µm | 5 / pk | 322.435 |
| PTFE Membrane filter, 0.2 µm | 5 / pk | 322.400 |
| Spare silicone pipette holder | 1 / pk | 1.435.01 |
| Spare nose cone in white | 3 / pk | 1.435.02 |
| Alternative nose cones in pink, violet and green (1 each) | 3 / pk | 1.435.03 |

For more information or for an overview of the complete Socorex product range, visit www.socorex.com