



Product Specification Sheet

TS/30-D

5 µl Red Culture Loops (soft), 20 in the sealed minigrip bag

Product Code	TS/30-D
Description	5 µl Red Culture Loops (soft), 20 in the sealed minigrip bag
Intended Use	Sample throughout. General laboratory use only
Components	5 µl Red Loops (soft) Minigrip Bag
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 20-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Laboratory use only
Pack Size	1000 in a box
Other	Quality control testing details on final product: Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 5 Date : 06.06.2016

ProductSpec_TS30-D