

(

# High Level Disinfectant Wipes

CHEMGENE HLD4I disinfectant wipes are manufactured with the latest innovation in nonwoven material, offering unrivalled strength, absorbency and low lint properties at a very cost effective price. The wipes are impregnated with a multi-active disinfectant solution that is alcohol free and tested to EN standards as effective against a range of bacteria, viruses, fungi, yeasts and mycobacteria.

# **BENEFITS OF USE ALCOHOL FREE**

### Compliant

- manufactured in accordance with Quality Management System accredited to ISO13485 and carry the CE mark of conformity
- Controlled dosage ensures effective disinfection with each wipe
- · Wipes are tested to a range of European Standards
- · Mechanical action enhances disinfection.

#### Convenient

- ·No requirement for pre-mixing, reducing wastage and risk of errors.
- · Canister can be stored on a work top or in a wall mounted bracket.

#### Cost effective

- · Saves valuable time and reduces the need for cloths and cleaning solutions.
- 100 wipes per canister means low cost per wipe.
- · Fully recyclable material and packaging



Class IIb registered







#### **AREAS OF USE**

These wipes are designed for use on high risk areas and on most surfaces, including glass, metal, ceramics, plastics and rubber. The tub can be stored in the area of use for easy access to quickly disinfect surfaces at risk of contamination.

#### **TECHNICAL INFORMATION**

CHEMGENE HLD<sub>4</sub>I Disinfectant wipes are developed and tested by qualified chemists and microbiologists. They are then subjected to stringent EN testing procedures to ensure efficacy against a wide range of Bacteria, Viruses, Fungi and Spores. Certificates are available on request.

Use biocides safely. Always read the label and product information before use.

#### PRODUCT SPECIFICATION

Wipe Material: 23gsm Polypropylene viscose

Quantity/Tub: 100 wipes
Outer Case Quantity: 10 tubs

Wipe Size:

100 Wipe Drum - 130mm x 180mm

Material Colours: White

#### CHEMGENE HLD | has been tested effective against:

#### **FUNGICIDAL ACTIVITY**

Trichophyton rubrum Aspergillus niger Candida albicans

#### VIRUCIDAL ACTIVITY

Hepatitis B HIV Herpes Simplex Norovirus

#### SPORICIDAL ACTIVITY

Bacillus subtilis Bacillus cereus

Coronavirus

## MYCOBACTERICIDAL

ACTIVITY

Mycobacterium terrae Mycobacterium avium

#### **BACTERICIDAL ACTIVITY**

Pseudomonas aeruginosa
Escherichia coli
Staphylococcus aureus
Enterococcus hirae
Legionella pneumophila
Salmonella typhimurium
Listeria monocytogenes
Methicillin-resistant
Staphylococcus aureus
Methicillin-resistant
Staphylococcus pseudointermedius

Further testing pending certification.

#### **QUANTITIES PER CASE**

100 wipe drum: 10 per case



#### **PURPOSE**

for the cleaning and disinfecting of non-surgically invasive devices such as transvaginal and transrectal probes / transducers and for hard surfaces of non-invasive, non-lumened medical devices sych as monitors, radiographic and diagnostic equipment.

Not to be used to disinfect surgically invasive medical devices





The details given in this leaflet are correct at the time of going to press. The company reserves the right to change and improve the products without prior notice.

Manufactured in the UK

For more information on our products or for technical support contact: East Point, High Street, Seal, Kent, TN15 0EG, UK

T 08452 223344 F 01732 763530

info@medi-mark.co.uk www.medi-mark.co.uk



