



## Product Specification Sheet

### **TS/5-1**

5ml Amies transport medium without charcoal in tube (blue cap) with blue cap cotton wool tipped(wood shaft) transport swab in a peel pouch

Product Code	TS/5-1
Description	5ml Amies transport medium without charcoal in tube (blue cap) with blue cap cotton wool tipped(wood shaft) transport swab in a peel pouch
Intended Use	This product is intended for bacteriological samples only. DO NOT USE THIS PRODUCT WHEN CHLAMYDIAL INFECTION IS EXPECTED
Components	Blue cap Polypropylene Tube Blue cap cotton/wood swab Non-toxic glue Label & film
Contents	De-ionised water Sodium chloride Disodium hydrogen phosphate Sodium thioglycollate Calcium chloride Magnesium chloride Potassium phosphate Agar
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Recovery after 24 hours will be unpredictable and target organisms will be lost. Fatty acid sensitive strains will give decreased recovery or be stressed in the presence of cotton
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 4 Date : 06.06.2016

ProductSpec\_TS5-1