



Product Specification Sheet

TS/5-22

5ml Stuart transport medium in tube with blue cap pure viscose tip polystyrene shaft swab in a peel pouch

Product Code	TS/5-22
Description	5ml Stuart transport medium in tube with blue cap pure viscose tip polystyrene shaft swab in a peel pouch
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Blue cap Polypropylene Tube Blue cap viscose/polystyrene swab Non-toxic glue Label & film Transport medium
Contents	De-ionised water Stuart Transport medium Agar
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterilityii) Organism Recovery

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 4 Date : 06.06.2016

ProductSpec_TS5-22