



Product Specification Sheet

TS/8-A

Sterile swab with a wooden shaft and cotton tip, single in the peel pouch.

Product Code	TS/8-A
Description	Sterile swab with a wooden shaft and cotton tip, single in the peel pouch.
Intended Use	Culture swab & general laboratory use
Components	Cotton tip/wooden shaft swab Peel pouch Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-46kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	Certain viruses and Micro-organisms are sensitive to cotton and maybe stressed or not even recover
Pack Size	1000 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 5 Date : 06.06.2016

ProductSpec_TS8-A