



## **Product Specification Sheet**

### **TS/15-T40**

Sterile individually wrapped blue Plastic Template 10x10 cm with press 'hold' tab

|                             |  |
|-----------------------------|--|
| Product Code                | TS/15-T40  |
| Description                 | Sterile individually wrapped blue Plastic Template 10x10 cm with press 'hold' tab                                    |
| Intended Use                | Defined area Sample Collection. General laboratory use only  |
| Components                  | Blue Plastic, easy open stomacher pouch.   |
| Contents                    | Not applicable.  |
| Method of Sterilisation     | Sterilised by gamma irradiation at dosage of 25-53.4kGy  |
| Shelf Life                  | 24 months from date of manufacture   |
| Storage Conditions          | 4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging |
| Precautions and Limitations | Template can only be used once or cross over contamination will occur  |
| Pack Size                   | 40 in a box  |
| Other                       | Quality control testing details on final product: i) Validation of irradiation sterility.                            |

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 5 Date: 24.04.2020  
ProductSpec\_TS15-T40