



## **Product Specification Sheet**

### **TS/15-YCS**

Sterile yellow sponge in easy open stomacher pouch, individually packed & dosed with 10 ml Peptone Saline Diluent (PSD)

Product Code	TS/15-YCS
Description	Sterile yellow sponge in easy open stomacher pouch, individually packed & dosed with 10 ml Peptone Saline Diluent (PSD)
Intended Use	Alternative material for sampling of meat carcasses for microbiological analysis
Components	Yellow Polyurethane Sponge, easy open stomacher pouch.
Contents	De-ionised Water – 1 Litre, Bacteriological Peptone – 1g, NaCl – 8.5 g.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	75 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism recovery (S.typhimurium)

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 2 Date: 24.04.2020  
ProductSpec\_TS15-YCS