



Product Specification Sheet

TS/165

Sterile 4000 ml Non-sterile Empty Sampling Container for 24 hour Urine Analysis

Product Code	TS/165
Description	Sterile 4000 ml Non-sterile Empty Sampling Container for 24 hour Urine Analysis
Intended Use	Multipurpose Use
Components	HDPE Bottle, HDPE Cap, Label.
Contents	Not applicable as supplied empty.
Method of Sterilisation	N/A - supplied non-sterile
Shelf Life	N/A
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	50 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 2 Date: 24.04.2020
ProductSpec_TS165