

**NEW !**



technoclone

## **CEVERON<sup>®</sup> 100 series**

First fully automated coagulation analyzer series for clotting, chromogenic, turbidimetric assays, thrombin generation and quenching technology (FXIII & ADAMTS13 activity)



*An explosion of possibilities !*

# The evolution of fully automated testing continues

**The first fully automated instrument series where you can combine classic haemostasis diagnostics with Thrombin Generation, Factor XIII activity\* and ADAMTS13 activity\***

Simultaneous measurement of five technologies.

Up to five technologies clotting, chromogenic, turbidimetric, thrombin generation and quenching can be run at the same time.

	Clotting	Chromogenic	Turbidimetric	Thrombin Generation	FXIII activity*	ADAMTS13 activity*
ceveron <sup>®</sup> c 1 0 0	✓	✓	✓			
ceveron <sup>®</sup> t 1 0 0	✓	✓	✓	✓		
ceveron <sup>®</sup> s 1 0 0	✓	✓	✓	✓	✓	✓



ADAMTS13  
FXIII activity  
Thrombin Generation  
chromogenic assays  
turbidimetric assays  
clotting assays